

STI Collection Kits and Substitutions 12.24.20 1400

Spectrum Health Contact: Laboratory Services *highlight denotes new content

Due to availability, we may have changed your order. There is currently a national shortage of collection devices. Please refer to the table below for order and collection information.

| Kit Name | Lawson # & Supply Availability | Replace supplies with if no inventory | Testing Performed | Collection Site/Method |
|---------------------------------------|--|---|--|--|
| Abbott Multi Collect Kit | #659762 NOT AVAILABLE | Depending on Specimen Collection Type: #450034 (Aptima unisex swab kit) or #80959 (Aptima multitest vaginal kit) or #80956 (Aptima first void urine) | Chlamydia Only PCR Urine (LAB260) Gonococcus Only PCR Urine (LAB2111089) Chlamydia + Gonococcus PCR Urine (LAB744) Chlamydia Only PCR Swab (LAB740) Gonococcus Only PCR Swab (LAB785) Chlamydia + Gonococcus PCR Swab (LAB2111533) | Urine, Endocervical, Vaginal, Male Urethral Cannot be used for children under 14yo Cannot be used for Rectal, Oral, Throat |
| Aptima Unisex Swab Kit | #450034 LIMITED Availability | #659762 Abbot Multi Collect Kit | Aptima Chlamydia + Gonococcus (LAB3522) Aptima Chlamydia NAAT (LAB3524) | Male Urethral, Endocervical Rectal, Oral, Throat sent to reference lab for testing |
| Aptima First Void Urine | #80956 LIMITED Availability | #659762 Abbot Multi Collect Kit | Aptima Gonococcus NAAT (LAB3523) Aptima Mycoplasma Genitalium (LAB1230560) Aptima Trichomonas NAAT (LAB3525) | Urine, Male or Female Children under 14yo sent to reference lab for testing |
| Aptima Multitest (Vaginal) Kit | #80959 Stocked | #659762 Abbot Multi Collect Kit | Aptima STI NAAT Panel (LAB1230566) | Vaginal Rectal, Oral, Throat sent to reference lab for testing |
| Thin Prep Container | Not available via Lawson. Order via Supply Order Form to Lab. Stocked | No alternative for replacement | Aptima Chlamydia + Gonococcus (LAB3522) Aptima Chlamydia NAAT (LAB3524) Aptima Gonococcus NAAT (LAB3523) Aptima Trichomonas NAAT (LAB3525) Pap Test (LAB1230097) HPV High Risk Screen (LAB263) HPV 16/18/45 Genotype (LAB3340) | Endocervical, Ectocervical |

Reminders:

- Always make sure vial is secured, leaky specimens may be rejected.
- **REQUIRED:** 2 patient identifiers on every specimen (full name and date of birth) and source.