

Office/Facility name		Date
Address	Ordered by Phone (	)

ID ORIGINAL WITH YOUR COURIER OR EMAIL TO: CHI	DUFFEIRE			
Y PHLEBOTOMY		QTY		
Tube gold 5 mL (SST)	(order as ea)		Advanced technology (ATL) requisition form	(100/pl
Tube red 6 mL (plain)	(order as ea)		Infectious disease (Microbiology) requisition form	(100/pl
Tube lavender 3 mL (EDTA)	(order as ea)		Allergen test requisition form	(100/pl
Tube light blue 2.7 mL (3.2% citrate)	(order as ea)		Surgical Pathology requisition	(100/pl
Tube dark green 6 mL (sod heparin)	(order as ea)		Additional test requisition (add on)	(order as
Tube light green 4.5 mL (lith heparin PST)	(order as ea)		Patient home collection instructions (English)	(order as
Tube royal blue 6 mL (trace element - Red stripe)	(order as ea)		Patient home collection instructions (Spanish)	(order as
Tube royal blue 6 mL MAYO (metal free EDTA - purple stripe)	(order as ea)		Orange Epic/Epic Care link lab order slips	(order as
Tube tan 3 mL (EDTA) or Tube royal blue 3 mL (EDTA)	(order as ea)		Purple Interface/E-Share slips	(order as
Tube pink 6 mL (EDTA)	(order as ea)	QTY	CYTOLOGY	
Tube gray 6 mL (sod fluoride)	(order as ea)		Pap test thin prep vial w/ cervical broom	(25/pl
Tube yellow 6 mL ACD (sol A)	(order as ea)		Pap test thin prep vial w/ cervical brush & spatula	(25/pl
Tube yellow 6 mL ACD (sol B)	(order as ea)		Rover Evalyn (self-collect) brushes with Thin Prep vials	(25/pl
Tube microtainer green (lith heparin)	(25/bag)		Cytolyt for non GYN small mouth	(order as
Tube microtainer gold (SST)	(25/bag)		Cytolyt for non GYN large mouth	(order as
Tube microtainer red (plain)	(25/bag)		Green top slide fixative	(order as
Tube microtainer lavender (EDTA)	(25/bag)		UroVysion kits	(order as
Tube QuantiFeron	(order as ea)	ОТУ	GLUCOLA	(0.00.00
Tube urine myoglobin (mayo)	(order as ea)		Glucola orange (50 gram)	(6/pl
Multi draw vacutainer needle 21 - 1½"	(48/box)		Glucola orange (100 gram)	(6/pl
Multi draw vacutainer needle 21 - 174  Multi draw vacutainer needle 22 - 11/4"	(48/box)		Glucola fruit punch (100 gram)	(6/p
Vacutainer holder	(=0/b0x)		Glucola lemon lime (50 gram)	(6/p
Clear aliquot tube	(25/bag)	OTV	SPECIMEN BAGS	(6/þ
Brown/Amber aliquot tube	(25/bag)	VII	,	(100/0
MAYO metal free aliquot tube	(25/bag)	_	White specimen bag (ambient)	(100/p
'		_	Green specimen bag (refrigerated)	(100/p
MAYO clear aliquot tube	(100/bag)	_	Blue specimen bag (frozen)	(100/p
TOXICOLOGY		_	Red specimen bag (stat)	(100/p
Blue top drug screen cups	(order as ea)		Purple specimen bag (priority)	(100/p
Drug screen cup seals (numbered)	(order as ea)		Biohazard specimen bag (large)	(order as
Drug screen seals (unnumbered)		QTY	MISCELLANEOUS	
Polymedco Fit Kit	(order as ea)		Breath test kit	(5/b
Y MICROBIOLOGY			Specimen log book	(order as
Standard tip Eswab	(order as ea)		Pediatric lead kit	(order as
Nasopharyngeal Eswab	(order as ea)		Pediatric lead kit (FILTER PAPER ONLY)	(order as
Blood culture bottle AFB	(order as ea)		72 hour stool container	(2/p
Blood culture bottle aerobic (green)	(order as ea)		Fetal fibronectin collection kit	(order as
Blood culture bottle anaerobic (orange)	(order as ea)		Panorama prenatal screen	(order as
Blood culture prep kit (Chloraprep)	(order as ea)		Horizon collection kit	(order as
Blood culture device	(order as ea)		Horizon/Panorama combo kit	(order as
Blood transfer device	(order as ea)		Cortisol saliva kit	(order as
Pinworm paddle	(order as ea)		Centrifuge (Conical) Tubes 15mL	(50/p
C & S, O & P stool kit	(order as ea)	OTY		(0.0)
STD multi-collection kit (Alinity)	(order as ea)	•	Container with formalin 20 mL	(order as
UTM viral transport tube	(order as ea)		Container with formalin 60 mL	(order as
Y URINE COLLECTION	(Order as ea)		Container with formalin 180 mL	(order as
	(	<b>—</b>		
Sterile orange top urine container - Max Limit 200	(order as ea)	_	Container with formalin 480 mL	(order as
24 hour urine container - Max Limit 10	(order as ea)	_	Container dry 180 mL	(order as
Urine hat - max limit 25	(order as ea)		Container dry 32 oz	(order as
Urine C & S kit (Gray tube with transfer straw)	(order as ea)		Container dry 83 oz	(order as
Pediatric urine collection kit (PUCK)	(order as ea)		Container dry 1 gallon/lids	(order as
Acetic acid (15 mL)	(6/pkg)		Flow Media/Hank solution	(order as
Acetic acid (25 mL)	(6/pkg)	L	Michels Solution	(order as
Boric acid 4 - 2.5 gram tablets	(order as ea)	QTY	LABELS	
	(order as ea)		Specimen label - Epic care link 1" core 300778	(r
Sodium carbonate 5 grams			Manual specimen labels	1)
	(order as ea)		Patient Registry 3x1 - 1" core 300770	(1
Sodium carbonate 5 grams	(order as ea) (order as ea)		1	,
Sodium carbonate 5 grams Sodium hydroxide 10 mL	,		Specimen label - Epic care link 3" core 300745	(1
Sodium carbonate 5 grams Sodium hydroxide 10 mL Urine yellow vacutainer 10 mL	(order as ea)	OTY		(1
Sodium carbonate 5 grams Sodium hydroxide 10 mL Urine yellow vacutainer 10 mL Urine transfer straw	(order as ea) (order as ea)	QTY	Specimen label - Epic care link 3" core 300745  OTHER	(r
Sodium carbonate 5 grams Sodium hydroxide 10 mL Urine yellow vacutainer 10 mL Urine transfer straw  REQUISITIONS, FORMS AND PATIENT EDUCAT	(order as ea) (order as ea)	QTY		(r
Sodium carbonate 5 grams Sodium hydroxide 10 mL Urine yellow vacutainer 10 mL Urine transfer straw   REQUISITIONS, FORMS AND PATIENT EDUCAT Laboratory supply requisition	(order as ea) (order as ea) (order as ea)	QTY		(r
Sodium carbonate 5 grams Sodium hydroxide 10 mL Urine yellow vacutainer 10 mL Urine transfer straw  REQUISITIONS, FORMS AND PATIENT EDUCAT	(order as ea) (order as ea)	QTY		(1

ALLOW 5 DAYS FOR DELIVERY OF SUPPLIES. Corewell Health Integrated Fleet Team (SHIFT) 616.486.5983