

Office/Facility name		Date	
Address	Ordered by	Phone ()	

ND ORIGINAL WITH YOUR COURIER OR EMAIL TO: CHLABSUPPLYRE PHLEBOTOMY		QTY	TY REQUISITIONS, FORMS AND PATIENT EDUCATION		
Tube gold 5 mL (SST)	(order as ea)	Ų.,	Advanced technology (ATL) requisition form	(100/pkg)	
<u> </u>		_	33 ((1 0/	
Tube red 6 mL (plain)	(order as ea)		Infectious disease (Microbiology) requisition form	(100/pkg)	
Tube lavender 3 mL (EDTA)	(order as ea)		Allergen test requisition form	(100/pkg)	
Tube light blue 2.7 mL (3.2% citrate)	(order as ea)		Surgical Pathology requisition	(100/pkg)	
Tube dark green 6 mL (sod heparin)	(order as ea)		Additional test requisition (add on)	(order as ea)	
Tube light green 4.5 mL (lith heparin PST)	(order as ea)		Patient home collection instructions (English)	(order as ea)	
Tube royal blue 6 mL (trace element - Red stripe)	(order as ea)		Patient home collection instructions (Spanish)	(order as ea)	
Tube royal blue 6 mL MAYO (metal free EDTA - purple stripe)	(order as ea)		Orange Epic/Epic Care link lab order slips	(order as ea)	
Tube tan 3 mL (EDTA) or Tube royal blue 3 mL (EDTA)	(order as ea)		Purple Interface/E-Share slips	(order as ea)	
Tube pink 6 mL (EDTA)	(order as ea)	ОТҮ	CYTOLOGY		
Tube gray 6 mL (sod fluoride)	(order as ea)		Pap test thin prep vial w/ cervical broom	(25/pkg)	
Tube yellow 6 mL ACD (sol A)	(order as ea)		Pap test thin prep vial w/ cervical brush & spatula	(25/pkg)	
Tube yellow 6 mL ACD (sol B)	(order as ea)	-	Rover Evalyn (self-collect) brushes with Thin Prep vials	(25/pkg)	
Tube microtainer green (lith heparin)	(25/bag)	_	Cytolyt for non GYN small mouth	(order as ea)	
9 1 1	. 0,		<u> </u>	, ,	
Tube microtainer gold (SST)	(25/bag)		Cytolyt for non GYN large mouth	(order as ea)	
Tube microtainer red (plain)	(25/bag)		Green top slide fixative	(order as ea)	
Tube microtainer lavender (EDTA)	(25/bag)		UroVysion kits	(order as ea)	
Tube QuantiFeron	(order as ea)	QTY	GLUCOLA		
Tube urine myoglobin (mayo)	(order as ea)		Glucola orange (50 gram)	(6/pkg)	
Multi draw vacutainer needle 21 - 11/4"	(48/box)		Glucola orange (100 gram)	(6/pkg)	
Multi draw vacutainer needle 22 - 11/4"	(48/box)		Glucola fruit punch (100 gram)	(6/pkg)	
Vacutainer holder	(each)		Glucola lemon lime (50 gram)	(6/pkg)	
Clear aliquot tube	(25/bag)	OTY	SPECIMEN BAGS	(7) 3/	
Brown/Amber aliquot tube	(25/bag)		White specimen bag (ambient)	(100/pkg)	
MAYO metal free aliquot tube	(25/bag)		Green specimen bag (refrigerated)	(100/pkg)	
MAYO clear aliquot tube	(100/bag)	_	Blue specimen bag (frozen)	(100/pkg)	
TOXICOLOGY	(100/bag)	_	Red specimen bag (stat)	(100/pkg)	
-	(-	1 3 ()		
Blue top drug screen cups	(order as ea)		Purple specimen bag (priority)	(100/pkg)	
Drug screen cup seals (numbered)	(order as ea)		Biohazard specimen bag (large)	(order as ea)	
Drug screen seals (unnumbered)		QTY	MISCELLANEOUS	<u> </u>	
Polymedco Fit Kit	(order as ea)		Breath test kit	(5/box)	
OTY MICROBIOLOGY			Specimen log book	(order as ea)	
Standard tip Eswab	(order as ea)		Pediatric lead kit	(order as ea)	
Nasopharyngeal Eswab	(order as ea)		Pediatric lead kit (FILTER PAPER ONLY)	(order as ea)	
Blood culture bottle AFB	(order as ea)		72 hour stool container	(2/pkg)	
Blood culture bottle aerobic (green)	(order as ea)		Fetal fibronectin collection kit	(order as ea)	
Blood culture bottle anaerobic (orange)	(order as ea)		Panorama prenatal screen	(order as ea)	
Blood culture prep kit (Chloraprep)	(order as ea)		Horizon collection kit	(order as ea)	
Blood culture device	(order as ea)		Horizon/Panorama combo kit	(order as ea)	
Blood transfer device	(order as ea)	_	Cortisol saliva kit	(order as ea)	
	, ,	_		, ,	
Pinworm paddle	(order as ea)		Centrifuge (Conical) Tubes 15mL	(50/pkg)	
C & S, O & P stool kit	(order as ea)	QIY	HISTOLOGY		
STD multi-collection kit (Alinity)	(order as ea)		Container with formalin 20 mL	(order as ea)	
UTM viral transport tube	(order as ea)		Container with formalin 60 mL	(order as ea)	
TY URINE COLLECTION			Container with formalin 180 mL	(order as ea)	
Storilo orango ton urino container Mary Limit 200	/ 1 \		Container with formalin 480 mL	(order as ea)	
Sterile orange top urine container - Max Limit 200	(order as ea)			, ,	
24 hour urine container - Max Limit 200	(order as ea)		Container dry 180 mL	(order as ea)	
			Container dry 180 mL Container dry 32 oz	(order as ea)	
24 hour urine container - Max Limit 10 Urine hat - max limit 25	(order as ea) (order as ea)		Container dry 32 oz	(order as ea)	
24 hour urine container - Max Limit 10 Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw)	(order as ea) (order as ea) (order as ea)		Container dry 32 oz Container dry 83 oz	(order as ea) (order as ea)	
24 hour urine container - Max Limit 10 Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK)	(order as ea) (order as ea) (order as ea) (order as ea)		Container dry 32 oz Container dry 83 oz Container dry 1 gallon/lids	(order as ea) (order as ea) (order as ea)	
24 hour urine container - Max Limit 10 Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK) Acetic acid (15 mL)	(order as ea) (order as ea) (order as ea) (order as ea) (6/pkg)		Container dry 32 oz Container dry 83 oz Container dry 1 gallon/lids Flow Media/Hank solution	(order as ea) (order as ea) (order as ea) (order as ea)	
24 hour urine container - Max Limit 10 Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK) Acetic acid (15 mL) Acetic acid (25 mL)	(order as ea) (order as ea) (order as ea) (order as ea) (6/pkg)	OEA	Container dry 32 oz Container dry 83 oz Container dry 1 gallon/lids Flow Media/Hank solution Michels Solution	(order as ea) (order as ea) (order as ea)	
24 hour urine container - Max Limit 10 Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK) Acetic acid (15 mL) Acetic acid (25 mL) Boric acid 4 - 2.5 gram tablets	(order as ea) (order as ea) (order as ea) (order as ea) (6/pkg) (order as ea)	QTY	Container dry 32 oz Container dry 83 oz Container dry 1 gallon/lids Flow Media/Hank solution Michels Solution LABELS	(order as ea) (order as ea) (order as ea) (order as ea) (order as ea)	
24 hour urine container - Max Limit 10 Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK) Acetic acid (15 mL) Acetic acid (25 mL) Boric acid 4 - 2.5 gram tablets Sodium carbonate 5 grams	(order as ea) (6/pkg) (6/pkg) (order as ea) (order as ea)	QTY	Container dry 32 oz Container dry 83 oz Container dry 1 gallon/lids Flow Media/Hank solution Michels Solution LABELS Specimen label - Epic care link 1" core 300778	(order as ea) (order as ea) (order as ea) (order as ea) (order as ea) (roll)	
24 hour urine container - Max Limit 10 Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK) Acetic acid (15 mL) Acetic acid (25 mL) Boric acid 4 - 2.5 gram tablets Sodium carbonate 5 grams Sodium hydroxide 10 mL	(order as ea) (order as ea) (order as ea) (order as ea) (ofder as ea) (6/pkg) (6/pkg) (order as ea) (order as ea) (order as ea)	QTY	Container dry 32 oz Container dry 83 oz Container dry 1 gallon/lids Flow Media/Hank solution Michels Solution LABELS Specimen label - Epic care link 1" core 300778 Manual specimen labels	(order as ea) (order as ea) (order as ea) (order as ea) (order as ea) (roll) (roll)	
24 hour urine container - Max Limit 10 Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK) Acetic acid (15 mL) Acetic acid (25 mL) Boric acid 4 - 2.5 gram tablets Sodium carbonate 5 grams Sodium hydroxide 10 mL Urine yellow vacutainer 10 mL	(order as ea) (order as ea) (order as ea) (order as ea) (of/pkg) (6/pkg) (order as ea)	QTY	Container dry 32 oz Container dry 83 oz Container dry 1 gallon/lids Flow Media/Hank solution Michels Solution LABELS Specimen label - Epic care link 1" core 300778 Manual specimen labels Patient Registry 3x1 - 1" core 300770	(order as ea) (order as ea) (order as ea) (order as ea) (order as ea) (roll)	
24 hour urine container - Max Limit 10 Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK) Acetic acid (15 mL) Acetic acid (25 mL) Boric acid 4 - 2.5 gram tablets Sodium carbonate 5 grams Sodium hydroxide 10 mL	(order as ea) (order as ea) (order as ea) (order as ea) (ofder as ea) (6/pkg) (6/pkg) (order as ea) (order as ea) (order as ea)		Container dry 32 oz Container dry 83 oz Container dry 1 gallon/lids Flow Media/Hank solution Michels Solution LABELS Specimen label - Epic care link 1" core 300778 Manual specimen labels Patient Registry 3x1 - 1" core 300770 Specimen label - Epic care link 3" core 300745	(order as ea) (order as ea) (order as ea) (order as ea) (order as ea) (roll) (roll)	
24 hour urine container - Max Limit 10 Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK) Acetic acid (15 mL) Acetic acid (25 mL) Boric acid 4 - 2.5 gram tablets Sodium carbonate 5 grams Sodium hydroxide 10 mL Urine yellow vacutainer 10 mL	(order as ea) (order as ea) (order as ea) (order as ea) (of/pkg) (6/pkg) (order as ea)	QTY	Container dry 32 oz Container dry 83 oz Container dry 1 gallon/lids Flow Media/Hank solution Michels Solution LABELS Specimen label - Epic care link 1" core 300778 Manual specimen labels Patient Registry 3x1 - 1" core 300770	(order as ea) (order as ea) (order as ea) (order as ea) (order as ea) (roll) (roll)	
24 hour urine container - Max Limit 10 Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK) Acetic acid (15 mL) Acetic acid (25 mL) Boric acid 4 - 2.5 gram tablets Sodium carbonate 5 grams Sodium hydroxide 10 mL Urine yellow vacutainer 10 mL Urine transfer straw	(order as ea) (order as ea) (order as ea) (order as ea) (of/pkg) (of/pkg) (order as ea)		Container dry 32 oz Container dry 83 oz Container dry 1 gallon/lids Flow Media/Hank solution Michels Solution LABELS Specimen label - Epic care link 1" core 300778 Manual specimen labels Patient Registry 3x1 - 1" core 300770 Specimen label - Epic care link 3" core 300745	(order as ea) (order as ea) (order as ea) (order as ea) (order as ea) (roll) (roll)	
24 hour urine container - Max Limit 10 Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK) Acetic acid (15 mL) Acetic acid (25 mL) Boric acid 4 - 2.5 gram tablets Sodium carbonate 5 grams Sodium hydroxide 10 mL Urine yellow vacutainer 10 mL Urine transfer straw TY REQUISITIONS, FORMS AND PATIENT EDUCAT	(order as ea) (order as ea) (order as ea) (order as ea) (of/pkg) (of/pkg) (order as ea)		Container dry 32 oz Container dry 83 oz Container dry 1 gallon/lids Flow Media/Hank solution Michels Solution LABELS Specimen label - Epic care link 1" core 300778 Manual specimen labels Patient Registry 3x1 - 1" core 300770 Specimen label - Epic care link 3" core 300745	(order as ea) (order as ea) (order as ea) (order as ea) (order as ea) (roll) (roll)	
24 hour urine container - Max Limit 10 Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK) Acetic acid (15 mL) Acetic acid (25 mL) Boric acid 4 - 2.5 gram tablets Sodium carbonate 5 grams Sodium hydroxide 10 mL Urine yellow vacutainer 10 mL Urine transfer straw TY REQUISITIONS, FORMS AND PATIENT EDUCATION Laboratory supply requisition	(order as ea) (order as ea) (order as ea) (order as ea) (of/pkg) (of/pkg) (order as ea)		Container dry 32 oz Container dry 83 oz Container dry 1 gallon/lids Flow Media/Hank solution Michels Solution LABELS Specimen label - Epic care link 1" core 300778 Manual specimen labels Patient Registry 3x1 - 1" core 300770 Specimen label - Epic care link 3" core 300745	(order as ea) (order as ea) (order as ea) (order as ea) (order as ea) (roll) (roll)	
24 hour urine container - Max Limit 10 Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK) Acetic acid (15 mL) Acetic acid (25 mL) Boric acid 4 - 2.5 gram tablets Sodium carbonate 5 grams Sodium hydroxide 10 mL Urine yellow vacutainer 10 mL Urine transfer straw TY REQUISITIONS, FORMS AND PATIENT EDUCAT	(order as ea) (order as ea) (order as ea) (order as ea) (of/pkg) (of/pkg) (order as ea)		Container dry 32 oz Container dry 83 oz Container dry 1 gallon/lids Flow Media/Hank solution Michels Solution LABELS Specimen label - Epic care link 1" core 300778 Manual specimen labels Patient Registry 3x1 - 1" core 300770 Specimen label - Epic care link 3" core 300745	(order as ea) (order as ea) (order as ea) (order as ea) (order as ea) (roll) (roll)	

ALLOW 5 DAYS FOR DELIVERY OF SUPPLIES. Corewell Health Integrated Fleet Team (SHIFT) 616.486.5983